



TacMed


İnş. San. Tic. Ltd. Şti.
Official Exclusive Distributor



Vathin[®] Single-use Flexible
Video Rhinolaryngoscope

VATHIN
Value Your Health

Vathin® Single-use Flexible Video Rhinolaryngoscope





Slim Diameter

A slim insertion tube allows easier approach to the nasal, pharyngeal and laryngeal cavities, making the scope of observation more comprehensive.



Field of View

The field can be thoroughly inspected with the scope's wide 110° field of view.



Bending Range

A $\pm 210^\circ$ wide angulation range allows a smoothly insertion into the upper airway and offers the flexibility to perform better during the procedure.

Specifications

Model	Color	Working channel ID(mm)	Insertion tube OD(mm)	Working length(mm)	Bending range	Field of view
RL-E1800	Grey	0	3.2	300	Up 0°~210° Down 0°~210°	110°± 15%
RL-E1E00	Green	2.2	4.9			



www.tacmedmedikal.com

info@tacmedmedikal.com

[Gaziler Caddesi Yenişehir Mah. No:341/G Konak-İZMİR](#)

[Tel : +90 232 449 00 47](tel:+902324490047) • [+90 533 386 26 91](tel:+905333862691)

