



**TacMed**  
İnş. San. Tic. Ltd. Şti.  
Official Exclusive Distributor



**TacMed**  
İnş. San. Tic. Ltd. Şti.  
Official Exclusive Distributor

[www.tacmedmedikal.com](http://www.tacmedmedikal.com)

**VATHIN** Vathin  
Medical  
Value Your Health

**VATHIN** Vathin  
Medical  
Value Your Health



Tel : +90 232 449 00 47 • +90 533 386 26 91

Gaziler Caddesi Yenişehir Mah. No:341/G Konak-İZMİR

[www.tacmedmedikal.com](http://www.tacmedmedikal.com)

[info@tacmedmedikal.com](mailto:info@tacmedmedikal.com)

[tacmedmedikal](https://www.instagram.com/tacmedmedikal) [tacmedmedikal](https://www.linkedin.com/company/tacmedmedikal) [tacmedmedikal](https://www.youtube.com/channel/UCtacmedmedikal)



Vathin® Single-use  
Flexible Ureteroscope/  
Cystoscope

Urology solution






# Vathin<sup>®</sup> Single-use Flexible Ureteroscope/ Cystoscope

Vathin<sup>®</sup> Single-use Flexible Ureteroscope / Cystoscope system is specially designed to meet the evolving needs of urology specialists, providing a comprehensive all-in-one-solution for healthcare professionals and offering a range of advantages that make it the ideal choice for clinical application. With its unique features, the system ensures smoother drainage and clearer picture quality, enabling urologists to perform procedures with utmost accuracy.



## Vathin<sup>®</sup> Single-use Flexible Ureteroscope



Model	Working length(mm)	Working Channel ID (mm/Fr)	Outer diameter of main hose OD (mm/Fr)	Bending section	Field of view	Bending lock	Suction function	Color
US-S170	700	1.2/3.6Fr	2.8/8.4Fr	U285°/D285°	120°±10%	Yes	Yes	Grey
Model	Working length(mm)	Working channel ID (mm/Fr)	Outer diameter of main hose OD (mm/Fr)	Distal end size (mm/Fr)	Bending section	LED number	Bending lock	Color
US-S2600	700	1.2/3.6Fr	2.5/7.5Fr	7.11Fr	U285°/D285°	2	Yes	Blue

## Main Features

>>



### Suction

Smoother drainage, clearer picture  
Easier to operate with rotatable suction tube.

### Bending lock

More convenient operation, more accurate positioning



### 90° Horizontal Rotation

Optimize the view without grip adjustment

### Riveted Snake bone

Providing flexibility & durability for better performance  
Bending capability  
Ureteroscope: 0~285° Cystoscope: 0~210°



## Vathin<sup>®</sup> Single-use Flexible Cystoscope



Model	Working length(mm)	Work Channel ID (mm/Fr)	Outer diameter of main hose OD (mm/Fr)	Bending Range	Field of view	Bending lock	Suction function
CY-S1B01	450	2 / 6Fr	4 / 12Fr	U210°	120°±10%	Yes	Yes
CY-S1B00				D210°		No	No